

---

## O<sub>2</sub> Gas Sensor

Order Code O2-BTA

### What's Included

- O<sub>2</sub> Gas Sensor
- 250 mL gas sampling bottle (Nalgene bottle with lid)

### Quick Start

1. Connect the sensor to the interface (LabQuest Mini, LabQuest 2, etc.).
2. Start the appropriate data-collection software (Logger *Pro*, Logger Lite, LabQuest App) if not already running, and choose New from File menu.

The software will identify the sensor and load a default data-collection setup. You are now ready to collect data.

If you are collecting data using a Chromebook™, mobile device such as iPad® or Android™ tablet, or a Vernier wireless sensor or interface, please see the following link for up-to-date connection information:

**[www.vernier.com/start/o2-bta](http://www.vernier.com/start/o2-bta)**

### Tips

- **Very important:** Do not place the sensor into any liquid. The sensor is intended only for measuring gaseous, not aqueous, O<sub>2</sub> concentration.

For full product documentation, see

**[www.vernier.com/manuals/o2-bta](http://www.vernier.com/manuals/o2-bta)**

**Note:** Vernier products are for educational use only.

---

## CO<sub>2</sub> Gas Sensor

Order Code CO2-BTA

### What is Included

- CO<sub>2</sub> Gas Sensor
- 250 mL gas sampling bottle (Nalgene bottle with lid)

### Quick Start

1. Set the range switch on the sensor.
2. Connect the sensor to the interface (LabQuest Mini, LabQuest 2, etc.).
3. Start the appropriate data-collection software (*Logger Pro*, *Logger Lite*, *LabQuest App*) if not already running, and choose New from File menu.

The software will identify the sensor and load a default data-collection setup. You are now ready to collect data.

If you are collecting data using a Chromebook™, mobile device such as iPad® or Android™ tablet, or a Vernier wireless sensor or interface, please see the following link for up-to-date connection information:

**[www.vernier.com/start/co2-bta](http://www.vernier.com/start/co2-bta)**

**Note:** Vernier products are for educational use only.

### Tips

- Choose the desired range on the sensor. Most experiments will use the Low range (0–10,000 ppm). The High range would be appropriate when measuring human respiration, as expired air is typically 40,000–60,000 ppm CO<sub>2</sub>.

- Allow the CO<sub>2</sub> Gas Sensor to warm up for about 90 seconds before collecting data. During warmup, the readings will be very low.
- Do not allow liquids to come in contact with the CO<sub>2</sub> Gas Sensor. The sensor is intended only for measuring gaseous, not aqueous, CO<sub>2</sub> concentration.

For full product documentation, see

**[www.vernier.com/manuals/co2-bta](http://www.vernier.com/manuals/co2-bta)**



**Measure. Analyze. Learn.™**  
**Vernier Software & Technology**

13979 S. W. Millikan Way • Beaverton, OR 97005-2886  
Toll Free (888) 837-6437 • (503) 277-2299 • FAX (503) 277-2440  
[info@vernier.com](mailto:info@vernier.com) • [www.vernier.com](http://www.vernier.com)

Rev. 2/4/16

Logger Pro, Logger Lite, Vernier LabQuest, Vernier LabQuest Mini, and other marks shown are our trademarks or registered trademarks in the United States.

All other marks not owned by us that appear herein are the property of their respective owners, who may or may not be affiliated with, connected to, or sponsored by us.



Printed on recycled paper

# Go! Link<sup>®</sup>

## (GO-LINK)



Go! Link is a single-channel USB interface used to connect many Vernier sensors directly to a Windows<sup>®</sup> or Macintosh<sup>®</sup> computer. Some Vernier sensors can also be used when Go! Link is connected to a Chromebook<sup>™</sup>. For a complete list of compatible sensors, visit [www.vernier.com/go](http://www.vernier.com/go)

## Getting Started

### Go! Link with a Computer

1. Install Logger Lite<sup>®</sup> or Logger Pro<sup>®</sup> software on your computer.
  - Logger Lite is free and can be downloaded at [www.vernier.com/loggerlite](http://www.vernier.com/loggerlite)
  - Logger Pro has many additional features and is available for purchase at [www.vernier.com/lp](http://www.vernier.com/lp)
2. Connect Go! Link to the computer's USB port and the sensor to the Go! Link.
3. Start the software. The sensor will be recognized and in most cases will display the live reading. Note: The thermometer shown below is only displayed when using a temperature probe.
4. Click Collect to begin collecting data, if desired, or open an experiment.


Double-click the column heading to edit

Double-click to change the graph title

Digital meter provides live readouts

Thermometer can show °C and °F

Grab axes to stretch



### Go! Link with a Chromebook

1. Install Graphical Analysis<sup>™</sup> (version 1.2 or newer, available in the Chrome web store in late spring 2015).
2. Connect Go! Link to the Chromebook's USB port and the sensor to the Go! Link.
3. Launch Graphical Analysis. The sensor will be recognized and in most cases will display the live reading.
4. Start data collection, if desired.

**NOTE:** Vernier products are designed for educational use. Our products are not designed nor are they recommended for any industrial, medical, or commercial process such as life support, patient diagnosis, control of a manufacturing process, or industrial testing of any kind.

## Specifications

USB specification	1.1
Maximum sampling rate	200 samples/second
Resolution	12 bit

## Related Products

### Go! Link Teacher Pack (order code: GL-TP)

Save money by purchasing a Teacher Pack of eight Go! Links. Includes free download of Logger Lite software. More information is available at [www.vernier.com/go](http://www.vernier.com/go)

### Go! Temp® (order code: GO-TEMP)

Go!Temp is an easy-to-use temperature probe that connects directly to the USB port of a computer, LabQuest, or Chromebook. Includes free download of Logger Lite software. More information is available at [www.vernier.com/go-temp](http://www.vernier.com/go-temp)

### Go! Motion® (order code: GO-MOT)

Go! Motion is a motion detector that connects directly to the USB port of a computer, LabQuest, or Chromebook. Includes free download of Logger Lite software. More information is available at [www.vernier.com/go-mot](http://www.vernier.com/go-mot)

## Warranty

Vernier warrants this product to be free from defects in materials and workmanship for a period of five years from the date of shipment to the customer. This warranty does not cover damage to the product caused by abuse or improper use.



**Measure. Analyze. Learn.™**  
**Vernier Software & Technology**

13979 S. W. Millikan Way • Beaverton, OR 97005-2886  
Toll Free (888) 837-6437 • (503) 277-2299 • FAX (503) 277-2440  
[info@vernier.com](mailto:info@vernier.com) • [www.vernier.com](http://www.vernier.com)

Rev. 1/12/15

Go! Link, Logger Pro, Logger Lite, Vernier LabQuest, Graphical Analysis, Go! Link, Go! Motion, and other marks shown are our trademarks or registered trademarks in the United States.

All other marks not owned by us that appear herein are the property of their respective owners, who may or may not be affiliated with, connected to, or sponsored by us.



Printed on recycled paper.



Tested To Comply  
With FCC Standards

FOR HOME OR OFFICE USE